




# ESR KIT


Easy Screw Removal KIT

**OSSTEM<sup>®</sup>**  
IMPLANT

- Excellent cutting performance of SR Drill
- Simple and quick removal in step 2
- Removal success rate over 90%




## Variety of Cases Possible for Prostheses Fails



Screw fracture




Guide ESR handle SR drill SR tip Torque handle



Abutment fracture




AR tip



Screw hex slip



AR tip (mini)



Internal thread crushing burden



Re-tap

**OSSTEM**<sup>®</sup>  
IMPLANT

**DEUTSCHE OSSTEM GmbH**

Mergenthalerallee 25, 65760 Eschborn, Germany

Tel: +49(0)6196 777 550 Fax: +49(0)6196 777 5529

E-mail: [info@osstem.de](mailto:info@osstem.de) [www.osstem.de](http://www.osstem.de)